MERIDIAN^{Classic} alpha

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 Achtung

 Trocken lagern

adeor medical AG
Martinshof 5 – 83626 Valley
Germany

CE ¹³⁰⁴
Tel. +49 8024 47717 0
Fax +49 8024 47717 99
eMail: office@adeor.com
www.adeor.com

EN

MERIDIAN^{Classic} alpha

Single-Use Perforators INSTRUCTIONS FOR USE

STERILE R **R_X Only** 

Failure to read and follow the warnings and the instructions below may result in serious patient injury or death, or malfunction.

 **Warnings:**

Although individual perforator models have been designed to be reliable, the user must be aware of the possibility that the perforator may not disengage.

Be careful to choose the right model for the intended use (see description below).

Never use the single-use perforator for more than one patient.

Always hold perforator **perpendicular (90°)** to the skull at point of penetration. Failure to observe (e.g. tipping/“stirring“) could cause product malfunction and consequent patient injury.

When cutting, allow the perforator to do the cutting with **controlled and constant pressure**. If the pressure is too low or too high, it may trigger a malfunction.

The perforator is intended for use on healthy bone. Use extreme caution when performing a skull perforation on patients with high intracranial pressure, attached dura and any other abnormalities such as unhealthy bone, and on infants, children and elderly patients.

Apply extreme caution when working on areas of thin bone, such as temporal bone.

Do not use the perforator in case of suboptimal anatomical conditions (e.g., unevenness).

Be especially careful at the end of the drilling process. The bones may be unstable and applying too much force may cause injuries. **Avoid leverage force to the perforator while it is placed in the cranium as this may lead to undesirable results, including but not limited to fracture of the hudson shaft.**

Do not use perforator in or near an existing perforation.

If the perforator stops during the drilling process or releases too early, it may not lock back in and a restart of the perforator may not be possible in the same hole.

Do not re-sterilize the perforator under any circumstances.

Do not use perforator if the packaging is damaged, opened and / or the expiration date has passed.

Intended Use:

The different Meridian™ single-use perforator models are self-locking perforators and are intended for perforating the cranium. Once the perforation is complete and the pressure is removed from the perforator, the device automatically releases if used correctly.

These instructions are intended for the following products:

- For skull bone thicknesses as thin as 3 mm: PER 7-11-S, PER 11-14-S, PER 11-14-SG
- For skull bone thicknesses as thin as 2 mm: PER C6-9-S, PER C8-11-S, PER C11-14-S

- For skull bone thicknesses as thin as 1 mm: PER 7-11-SP, PER 11-14-SP

Inspection:

To be performed prior to each use:

- Check the expiry date printed on the label.
- Conduct a visual inspection of the packaging and perforator prior to opening.
- Open the packaging carefully.
- Inspect the perforator for any visible damage. Do not use the perforator if it shows any signs of damage.
- Perform a function test, as described below.

Function test:

This test is to be performed prior to each use. Do not use the perforator if any of the previous steps are not possible.

Use gloves while performing the following tests Only use perforators that run smoothly and lock/release automatically. The function test should be performed as follows:

- If the perforator is not pressed together, the blade and the shaft must be rotated against each other.

- It must be possible to press the shaft of the perforator and the blade against each other (see photo).

Operating Procedure:

- The shaft of the perforator fits into all conventional Hudson fittings.
- Insert the perforator securely into the drill system. Check that the connection between the perforator and the Hudson fitting is secure.
- Check that the settings on the control unit for the drill system have been set correctly. The authorized speed range is **800-1200 rpm**. The drill system must only rotate in a **clockwise** direction. **Counterclockwise drilling with the perforator may lead to undesirable results, including but not limited to overheating and disassembly.**
- To start trepanation, the perforator must be placed absolutely perpendicular (90°) to the calotte of the skull and must not be moved back and forth.
- Twist the perforator driver manually until the mechanism initially locks.
- Start the trepanation. Only consistent and controlled pressure may be applied during this process. Once sufficient perforation has been achieved, the perforator blade will stop. The surgeon must not exert any further pressure, as this can result in injury of the dura.

Repeat function test before each perforation.

The surgeon must be aware of the possibility of the skull bones being completely perforated before the release mechanism stops the rotation movement of the blade. If this happens, the surgeon must stop drilling immediately to avoid injuring the dura. **As soon as the auto-stop mechanism of the perforator engages and the blades stop rotating, the motor must be switched off. Continued operation may lead to undesirable results, including but not limited to overheating and disassembly.**

 **Re-Sterilization:**

The product in question is **only intended for single use**. It must not be re-sterilized. Always use an aseptic approach. Re-sterilization may cause product components to malfunction, which could negatively impact the function of the product.

 **Storage & Handling:** in dry area.

Technical Data:

Rotation speed: 800-1200 rpm
Material: medical stainless steel
Sterility: Gamma-sterilized

Disposal:

Dispose of used and contaminated perforators immediately according to the standard hospital procedure for disposing of sharp, contaminated surgical instruments.

Service:

For service and support, please contact your local adeor® distributor. Unauthorized modifications carried out on medical equipment will result in loss of warranty rights.

Symbols used:

 Manufacturer

 Refer to user manual

R_X Only Caution: Federal law restricts this device to use by or at the order of a physician.

 Expiration Date

 Single-Use only

 Do not use if damaged or opened

STERILE R Sterilized using gamma irradiation

LOT Lot number

REF Catalogue number

 Caution

 Keep dry

adeor medical AG
Martinshof 5
83626 Valley
Germany

CE ¹³⁰⁴
Tel. +49 8024 47717 0
Fax +49 8024 47717 99
eMail: office@adeor.com
www.adeor.com

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