CRANIOBRIDGE

Screwdriver (handle and blade)
Schraubendreher (Griff und Klinge)
Cacciavite (impugnatura e lama)

- EN** Instructions for use
- DE** Gebrauchsanweisung
- IT** Istruzioni per l'uso

Table of Contents / Inhaltsverzeichnis / Indice

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SCREWDRIVER (handle and blade)

General information

Please read the following information carefully. Improper handling, maintenance, and inappropriate use can cause premature wear or risks for patients and user. This Instruction for Use does not contain all the information required to use a surgical instrument.

adeor medical AG products may only be used with original adeor medical AG equipment. The operating physician / surgeon must ensure that the equipment is in faultless working order prior to using it. Instruments have a limited service life; this lifespan can be reduced even further through the lack of care in handling or in case of inadequate protection. Any instrument that does not work properly (e.g. due to wear and tear, mishandling, improper maintenance) must be disposed of.

Intended use

The aim of using the instruments is the osteosynthesis in different anatomical regions. Thereby it is possible to fix bone fragments minimally invasive, selective and without compression effect. The instruments are intended for the insertion of the implants and are intended to support the recovery.

Sorting out damaged instruments

Check the instruments for damage and wear. Also check that the screws can be securely attached to the blade. Damaged instruments must not be used under any circumstances. Non-observance can lead to instrument breakage and endanger the patient and the user.

Combination with other products

If instruments are reassembled after disassembly, individual parts may not be replaced by parts from other manufacturers. Parts from other manufacturers must not be used even if parts are interchangeable (e.g. various applications) due to the construction of the product.

Precautionary measures before use

Upon delivery, the instruments are clean but not sterile. It is the responsibility of the hospital to decontaminate, clean and sterilise the instruments before use through validated processes. The following recommendations do not replace the current health regulations

(standards, guidelines, national recommendations, official documents, experience, etc.)

Care and maintenance are mandatory and necessary for preservation and functionality.

Preparation of the instruments must be carried out according to these instructions and is based on a validated process.

Hygiene and Preparation

Prior to any surgery, the packaging should be removed to perform a visual inspection and to ensure that all instruments are in proper condition. Instruments consisting of removable parts should be dismantled before decontamination and cleaning.

- ▶ Factory-new instruments and instruments coming back from repairs / returns should be prepared and processed in the same way as used instruments before their first use. Protective transportation packaging, protective caps, etc. are not suitable for sterilisation.
- ▶ Only approved agents (e.g. RKI, DGHM/VHA, FDA, etc.) may be used.
- ▶ Both alkaline and enzymatic cleaning agents are applicable.
- ▶ Water quality should be according to DIN EN 285 Annex B. The steriliser should be according to DIN EN 285 or DIN EN 13060.
- ▶ Cleaning and disinfection equipment should be according to DIN EN ISO 15883 Parts 1 and 2.
- ▶ Only appropriate equipment and product-specific validated procedures may be used for cleaning / disinfection / sterilisation.
- ▶ The instructions and recommendations of the manufacturer must be observed.
- ▶ Furthermore, the applicable national statutory / hygienic regulations must be observed. This applies especially to the various specifications regarding effective prion inactivation.

Processing

The instruments must be properly disassembled and machine cleaned before first use and after each use. If processing is to be performed manually, the individual parts of the instruments must be placed in a cleaning-active disinfectant. In this respect, please follow the instructions of the disinfectant manufacturer. All surfaces - including the internal cavities, lumens and openings - should come into contact with the solution. adeor recommends the following procedure for the cleaning and processing of the instruments:

Decontamination

Decontamination and cleaning are obligatory prior to sterilisation. The aim of the decontamination is to reduce the bacterial population, to make cleaning

easier at a later stage, to protect personnel and to avoid any contamination of the environment.

Immediately after use, instruments should be disassembled and immersed in distilled water or covered on a tray with a damp cloth.

To prevent the corrosive effect of salines, they should be removed immediately after surgery. The instruments must be decontaminated as soon as possible after use. For this purpose, the instruments should be completely immersed in an aldehyde-free, decontaminating, pH-neutral solution (proteolytic enzyme solution). The decontaminating solution must have a temperature of less than 30 °C to prevent the adhesion of bacteria. Soak times and rinsing technique must be followed according to the precautions for use of the decontaminant.

Please note that certain agents may cause discolouration and/or corrosion of the instruments. Examples of such agents include bleach, aqueous hypochlorite solution and sodium chloride. It is strongly recommended to soak the instruments in an ultrasonic cleaner. The instruments must then be rinsed out carefully. The decontamination solution must be replaced after each use. The use of some products may leave phosphates on the instruments after drying. To avoid corrosion of all kinds: do not leave the instruments in the decontamination bath too long and dry the instruments thoroughly.

Mechanical Processing

Important information on mechanical processing

- ▶ Remove coarse contaminants from the instruments immediately after use.
- ▶ Do not use fixing agents or hot water (> 40 °C) as this will cause residues to stick and may affect cleaning performance.
- ▶ The instruments must be disassembled or opened for processing.
- ▶ Wire brushes, scouring pads or other objects which could damage the instruments may not be used.
- ▶ Cleaning agents used must have a neutral pH to avoid corrosion and blunting of the instruments, and to prevent the risk of skin irritation of the user.
- ▶ The precautions for use as well as the recommendations of the washer-disinfector (WD) supplier should be observed regarding cleaning agents and washing cycles.
- ▶ The cleaning cycle must undergo a final rinse at a temperature set between 70 °C and 95 °C to facilitate drying. The instruments should be dried properly to prevent new contamination.
- ▶ After cleaning, each instrument should be inspected to see if it handles well and operates properly.

Pre-cleaning

- ▶ Place the instruments in cold water for at least 5 min. If possible, clean them under cold water

with a soft brush until no more residues are visible. In the case of cavities, bores and threads, rinse with a water gun for at least 10 sec (pulsed procedure).

- ▶ Place instruments in an ultrasonic bath for 15 min at 40 °C with a 0.5 % alkaline cleaning agent and treat with ultrasound.
- ▶ Remove the instruments and rinse with cold water.

Processing

Place the disassembled instruments on the inset of the cleaning tray. If they cannot be plugged onto the cleaning tray place them in a sieve tray.

- ▶ Pre-rinse with cold water for 1 min.
- ▶ Drain.
- ▶ Pre-rinse with cold water for 3 min.
- ▶ Drain.
- ▶ Wash for 5 min at 55 °C with a 0.5 % alkaline cleaning agent.
- ▶ Drain.
- ▶ Neutralisation for 3 min with warm tap water (> 40 °C) and a neutralising agent.
- ▶ Drain.
- ▶ Intermediate rinsing for 2 min with warm tap water (> 40 °C).
- ▶ Drain.

Disinfection

Perform mechanical thermal disinfection considering the national requirements regarding the A₀ value (see ISO 15883).

Disinfection parameters	
Temperature	93 °C (A ₀ 3000)
Soaking time	300 s (5 min)
Drying temperature	100 °C
Drying time	20 min

Sterilization

Products should be sterilised using a fractionated pre-vacuum process (according to ISO 13060 / ISO 17665) considering the respective national requirements.

- ▶ 3 pre-vacuum phases with at least 60 millibar pressure heating
- ▶ Sterilization temperature of at least 134 °C and a maximum of 137 °C
- ▶ Shortest hold time: 5 min
- ▶ Shortest drying time: 10 min

It is the responsibility of the user to ensure that the processing process, including resources, materials and personnel, is adequate for achieving the required results. State of the art technology and national laws require compliance with validated processes.

Pre-vacuum sterilization parameters (three phases)	
Pressure	> 60 mbar
Temperature	134 °C – 137 °C
Holding time	> 5 min
Drying time	> 10 min

If there is residual moisture in the sterilization tray after the sterilization cycle, proceed as follows:

- ▶ Do not open the box directly after use.
- ▶ Increase the drying time.

The recommendations for sterilization are only applicable as a guideline. Under no circumstances can the manufacturer be held responsible for the sterility of devices sterilized in the hospital.

- ▶ It is the responsibility of the hospital to decontaminate, clean and sterilize the instruments (complete box or individual instruments) in conformity with the recommendations of the manufacturer adeor.
- ▶ Any lost or damaged instruments should be reported to the responsible person, to the supplier, or directly to the manufacturer.

Sort out damaged instruments.



Packaging

- ▶ Packaging may be used pursuant to DIN EN 868.
- ▶ The packaging must be selected in such a way that the instruments fit into the packaging.
- ▶ Use a sterilization indicator for the packaging and note the sterilization and expiry date on the packaging.

Storage

- ▶ Sterile instruments should be stored dry, dust-free, without external impact, without major temperature fluctuations and not in the immediate vicinity of corrosive substances.
- ▶ It is advisable to store the instruments in trays, containers and cabinets.

Symbols used / explanation

	Consult instructions for use
	Manufacturer
	Date of manufacture
	Reference number
	Serial number

	Batch number
	Non-sterile
	Only use if packaging is undamaged
	WARNING: Observe warnings
	CE mark
	Medical device

Warranty

The products of adeor medical AG are developed and manufactured with greatest care.

Tests and inspections by highly qualified personnel guarantee safe and perfect functioning. Please note that warranty claims are only valid if all instructions in the enclosed instructions for use are followed.

adeor medical AG is not liable for damage caused by improper handling or repair by third parties not authorized by adeor.

The use of the instruments is the responsibility of the user. Any liability on the part of adeor medical AG is excluded for any resulting damage.

Warranty claims must be submitted to adeor or the service partner, enclosing the proof of purchase.

The provision of warranty service does not extend the warranty period or any other guarantee period.

adeor medical AG
 Martinshof 5
 83626 Valley
 Germany

Phone +49 8024 47717 0
 Fax +49 8024 47717 99
 Email: office@adeor.com
<http://www.adeor.de>

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I dispositivi di adeor medical AG sono sviluppati e realizzati con la massima cura.